

TRIAL TRANSFER CATHETER

Description:

CODE	Description
T-EB-ETC18	Trial Transfer Catheter for Embryo Trans (18 cm)
T-EB-ETC23	Trial Transfer Catheter for Embryo Trans (23 cm)
T-ET-ETC18	Trial Transfer Catheter for Echo Trans (18 cm)
T-ET-ETC23	Trial Transfer Catheter for Echo Trans (23 cm)
T-TEB-ETC18	Trial Transfer Catheter for Twinkle Embryo Trans (18 cm)
T-TEB-ETC23	Trial Transfer Catheter for Twinkle Embryo Trans (23 cm)

Indication for use:

Allwin trial transfer catheter is a closed ended catheter device provided sterile for single-use only. Allwin trial transfer catheters are only to be used to assess the passage through the cervix prior to embryo transfer to determine whether the cervix is passable for Allwin embryo replacement catheter

Contraindications:

This catheter should not be used on a patient with an active vaginal or intrauterine infection, a intra – fallopian procedure, a sexually transmitted disease, uterine perforation, a recent pregnancy (or is currently pregnant), or confirmed or suspected intrauterine device.

Complications:

- Urinary tract Infection
- Pelvic inflammatory disease
- Uterine infection
- Bleeding
- Endometrial / Endocervical damage
- Endometrial lesions

Caution:

- Federal (USA) law restricts this device to sale by or on the order of a physician.
- The visibility of catheter may be compromised with use of low-resolution ultrasound equipment.
- The patients should have a full bladder for performing transabdominal ultrasound.

Warning / Precaution:

- Do not use if package is opened or damaged.
- To reduce the risk of perforation, if any resistance is felt while inserting the stylet, soft obturator, or catheter, do not force the device against the resistance.
- Confirm the specification (length) of the catheter prior to use.
- This device is supplied sterile and intended for single use only.
- Do not re-sterilize.
- Only experienced practitioners should use this device.
- Infection may occur due to bacterial contamination of the device during vaginal manipulation, and result in urinary tract infection (UTI), pelvic inflammatory disease (PID), or uterine infection.
- Recommendations to minimize infection include the use of only embryo compatible materials, flushing the catheter (any other accessories used) with sterile, compatible culture media, and closely adhering to sterile techniques.
- Bleeding may occur as a result of the trauma due to insertion of the catheter through the cervix and has been reported to be associated with a lower pregnancy rate. A simple and atraumatic transfer method has been noted to offer the best condition for success.
- Dispose in accordance with recognized medical practice and under observance of applicable law and regulation.

- Storage temperature between 5°C to 30°C.

Instruction for Use:

1. Place patient in a lithotomy.
2. Insert vaginal speculum to expose cervix and clean with a cotton swab moistened with normal saline or medium.
3. Shape the outer catheter to complement the patient's anatomy and pass the catheter till it reaches the internal orifices.
4. Advance the catheter so that the inner catheter passes through the external and internal os, into the mid-uterine cavity. (It may be necessary to twist the catheter as it negotiates the cervical canal)
5. In case of difficulty in negotiating the cervical canal, use appropriate stylets as under

LENGTH	FOR USE WITH	CODE
18 cm	T-EB-ETC18, T-ET-ETC18, T-TEB-ETC 18	STY18
23 cm	T-EB-ETC23, T-ET-ETC23, T-TEB-ETC 23	STY23

6. Place the tip of catheter approximately 1cm from the fundus. Assess the passage of the catheter in preparation for embryo transfer.
7. Remove and dispose of the catheter in accordance with local medical hazardous waste practices.

Limited Express Warranty:

The limited express warranty as set forth herein is exclusive and in lieu of all warranties of merchantability and fitness for use, remedies, obligations and liabilities for consequential damages.

The products are being sold only for the purpose described herein and such limited express warranty runs only to the original user. In no event shall ALLWIN be liable for any breach of warranty in any amount exceeding the purchase price of the product. ALLWIN reserves the right to make design changes to products without liability to incorporate said changes in ALLWIN products previously designed or sold.

Note

Bacterial endotoxin (LAL - limulus amoebocyte lysate) assay ≤ 20 EU/device. Test conducted to ensure efficiency and safety of the device(s) for its intended use.



Read Instruction



For Single Use Only



Do not Re-Sterilize



Use By



Sterile unless package is opened or damaged



Caution : Federal (USA) law restricts device to sale by or on the order of a physician



Sterilized using ethylene oxide



Keep out of sunlight



Keep Dry



Not made with natural rubber latex



allwin Medical Devices

Allwin Medical Devices

3305 E. Miraloma Ave., Suite 176 Anaheim, CA 92806 (USA)

Tel. : +1 714-572-1709 | Fax : +1 714-572-1732

Email : info@allwinmedical.com | www.allwinmedical.com