

## Stone Basket & Extractors

### Description:

Trapper No Tip Nitinol Stone Basket  
Grabit Stainless Steel Stone Basket – Flatwire  
Grabit Stainless Steel Stone Basket – Helical  
Nitiflex Nitinol Stone Basket – Helical  
Nitiflex Nitinol Stone Basket – Flatwire  
Maestro Nitinol Tipless Stone Basket  
Elite Stone Retriever

### Configuration:

See brochure

### Intended use:

These devices have been designed for endoscopically manipulating, dislodging, engaging and removing renal, ureteral or biliary stones.

### Contraindications:

Not suited for use as a lithotripter basket for manual lithotripsy.

### Complications:

Some possible complications may include, but are not limited to, the following:

- Perforation
- Evulsion
- Edema
- Entrapment
- Hemorrhage
- Inability to disengage from irretrievable object

### Caution:

Federal (USA) law restricts this device to sale by or on the order of a physician.

### Warning/ Precaution:

- Manipulate the handle prior to use to ensure smooth functioning of the device. Avoid contact with electric instruments.
- Do not move the stone basket against sharp object.
- Check compatibility with the endoscope.
- Only physicians trained in stone manipulation should perform this procedure. A variety of techniques can be employed; however the physician should employ the technique with which he/she is most comfortable.
- For single use only.
- Attempt to resterilize a single use device may compromise its structural integrity and/or lead to device failure which may result in patient injury or illness.
- Do not use if package is opened or damaged.
- Storage temperature between 5°C to 30°C.

### Potential Biohazard:

Dispose in accordance with recognized medical practice and under observance of applicable law and regulation.

- Manipulate the handle to open the instrument slowly. Once the stone has been captured, partially close the instrument to secure the stone for removal.
- If unsuccessful, further attempts may be made by repeating previously described procedure.
- Simultaneously, withdraw the instrument and the endoscopes from the area of interest.

### Do not reuse, reprocess & resterilize:

Reuse, Reprocessing & resterilization may compromise the structural integrity, can also create risk of contamination & may not give desired result or create complications, infections which may result in injury, illness or death.

### Limited Express Warranty:

The limited express warranty as set forth herein is exclusive and in lieu of all warranties of merchantability and fitness for use, remedies, obligations and liabilities for consequential damages.

The products are being sold only for the purpose described herein and such limited express warranty runs only to the original user. In no event shall ALLWIN be liable for any breach of warranty in any amount exceeding the purchase price of the product. ALLWIN reserves the right to make design changes to products without liability to incorporate said changes in ALLWIN products previously designed or sold.

### Instruction for Use:

- Retract the basket / forceps into the catheter shaft.
- Insert the catheter shaft into the working channel of the endoscope.
- With the tip still retracted, advance the instrument towards the stone to be removed, which is fluoroscopically / endoscopically verified.



Catalogue  
Number



Batch  
Code



Date of  
Manufacture



Use  
By



Do not  
re-use



Do not  
Re-Sterilize



Caution: Federal Law restricts this device to sale by or on the order of a physician or a practitioner trained in its use.



Consult  
Instruction  
for use



Do not use if  
Packing is damaged



Caution



Keep out  
of sunlight



Keep Dry



Sterilized using  
ethylene oxide



CMC Medical Devices & Drugs S. L.  
C/Horacio Lengo N° 18, CP 29006, Malaga, Spain  
Tel.: +34 951 214 054  
E-mail : info@cmcmedicaldevices.com



**Allwin Medical Devices, Inc.**  
3305 E. Miraloma Ave., Suite 176 Anaheim, CA 92806 (USA)  
Tel. : +1 714-572-1709  
E-mail : info@allwinmedical.com | www.allwinmedical.com