

Tumor Stent

Description:

Tumor Stent

Configuration:

See brochure

Intended use:

Tumor Stents are used for temporary as well as for palliative treatments of internal urinary drainage when tumors are present in the ureteral tract

Tumor Stent is specially designed to optimize drainage, even in cases where the ureter is extremely compressed.

Contraindications:

- Untreated infected urinary congestion of the kidneys (pyonephrosis).
- Vascular abnormalities of the renal pelvis outlet.
- Acute urethral trauma.

Complications:

The potential complications associated with the use of the stent which includes, but are not limited to:

- Encrustation
- Hemorrhage
- Infection
- Stent Migration
- Perforation

Caution:

Federal (USA) law restricts this device to sale by or on the order of a physician.

Warning / Precaution:

- Only experience surgeons should place this device.
- Periodic evaluation is recommended. Not intended for permanent use.
- Check compatibility of guidewire before use.
- Do not use if package is opened or damaged.
- Dispose in accordance with recognized medical practice and under observance of applicable law and regulation.
- For single use only.
- Attempt to resterilize a single use device may compromise its structural integrity and/or lead to device failure which may result in patient injury or illness.
- Storage temperature between 5°C to 30°C.

Instruction for Use:

The Super Strong Tumor stent is advanced over a guidewire into the previously dilated ureter either retrograde or antegrade depending on the physician's choice.

Stent Removal:

1. Either using a suture or cystoscopically, retrieve the stent by gentle pulling.
2. If any resistance is encountered during removal, fluoroscopically determine stent position and cause of resistance.

Storage:

The device should be stored in a clean, dry area. Do not expose to organic solvents, ionizing radiation or ultraviolet light. The devices must be used prior to the sterilization expiration date on package.

How Supplied:

The Tumor Stent is provided sterile by ethylene oxide and is intended for single use only. Do not resterilize, reuse or reprocess.

Do not reuse, reprocess & resterilize:

Reuse, Reprocessing & resterilization may compromise the structural integrity, can also create risk of contamination & may not give desired result or create complications, infections which may result in injury, illness or death.

Limited Express Warranty:

The limited express warranty as set forth herein is exclusive and in lieu of all warranties of merchantability and fitness for use, remedies, obligations and liabilities for consequential damages. The products are being sold only for the purpose described herein and such limited express warranty runs only to the original user. In no event shall ALLWIN be liable for any breach of warranty in any amount exceeding the purchase price of the product. ALLWIN reserves the right to make design changes to products without liability to incorporate said changes in ALLWIN products previously designed or sold.



Catalogue
Number



Batch
Code



Date of
Manufacture



Use
By



Do not
re - use



Do not
Re-Sterilize

Rx only

Caution : Federal (USA) law
restricts device to sale by or
on the order of a physician.



Consult
Instruction
for use



Do not use
if Packing
is damaged



Caution



Keep out
of sunlight



Keep
Dry

STERILE EO

Sterilized using
ethylene oxide



30°C
Temperature
Limit



allwin
Medical Devices

Allwin Medical Devices, Inc.

3305 E. Miraloma Ave., Suite 176 Anaheim, CA 92806 (USA)

Tel. : +1 714-572-1709

E-mail : info@allwinmedical.com | www.allwinmedical.com