

Product Catalog

WOMEN'S HEALTH

🕲 www.allwinmedical.com



ITRODUCTION

OVUM PICKUP NEEDLES

ACE Single Lumen Ovum Pickup Needle	3
ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line	
ACE-M Manual Single Lumen Ovum Pickup Needle	
OVUMPICK Double Lumen Ovum Pickup Needle	

EMBRYO TRANSFER CATHETERS

EMBRYO TRANS Embryo Transfer Catheter	7
ECHO TRANS Embryo Transfer Catheter	
TWINKLE EMBRYO TRANS Embryo Transfer Catheter	
Trial Transfer Catheter	
BULB TRANS ULTRA Embryo Transfer Catheter - Bulb Tip	11
BULB TRANS STAR Embryo Transfer Catheter - Bulb Tip.	
PRO EMBRYO TRANS Embryo Transfer Catheter	
PRO ECHO TRANS Embryo Transfer Catheter	14
TWINKLE PRO EMBRYO TRANS Embryo Transfer Catheter	

INTRA UTERINE INSEMINATION CATHETERS

INSEMINA Intra Uterine Insemination Catheter	16
IUI CURVED Intra Uterine Insemination Catheter Curved	17
TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet	18
FLEXI Intra Uterine Insemination Catheter	
STUMPY Intra Uterine Insemination Cannula	20

FALLOPIAN TUBE RECANALIZATION CATHETERS

HYSTEROSALPINGOGRAPHY CATHETER

ENDOMETRIAL BIOPSY SAMPLING SYSTEM

THE SAMPLER Endometrial Biopsy Catheter	23
PRO-VAC Endometrial Biopsy Sampler	23

VITRIFICATION STRAWS

SURESTRIP Vitrification Straws (Open System Vitrification)	
SURELOCK Vitrification Straws	
SUREHOOK Vitrification Straws	

OTHER PRODUCTS



Product Warranty

Allwin warrants that at the time of manufacturing, the products are prepared and tested in accordance with highest standards. Due to biological differences in individuals, no product is 100% effective under all circumstances for all patients. Further, because we have no control of the condition under which the product is used, diagnosis of patient, the method of use or administration & handling and storage of the product after it leaves our possession, we do not warrant any consequence which may arise pertaining to the final use of the device. The foregoing warranty is exclusive and in lieu of all other warranties written, oral or implied. No representative of Allwin may change any of the forgoing and the buyer accepts the product subject to all terms hereof.

Design Changes

Continuous improvements is the chosen path at **Allwin** and we reserve the right to change or discontinue any product without notice.

Return Goods Policy

Allwin will accept returned goods only after issuing a Return Goods Authorisation (RGA). The return request must provide a reason for return, invoice number and order number. Credit will be issued once the request is evaluated as per company rules. We will accept a return of the stock item within 30 days of the original invoice date. We will not accept the return of a nonstock, opened and damaged products. Returned stock items may be subject to a restocking fee.

Product testing

Allwin products are endotoxin tested with a limulus amebocyte lysate (LAL) assay unless otherwise noted. Our Ovum Pickup needles and Embryo Transfer catheters are batch tested for toxicity with a one-cell mouse embryo assay (MEA) and our Intra Uterine Insemination catheters are Human Sperm Survival Assay (HSSA) tested.

Storage of Products

Allwin recommends that products be stored in a cool and dry place. Plastic material used may degrade following extended periods of exposure to light and temperature.




ACE Single Lumen Ovum Pickup Needle



The ACE Single Lumen Ovum Pickup Needle is used for laparoscopic or ultrasound-guided transvaginal retrieval of oocytes from ovarian follicles.

- Highly echogenic tip
- Silicone bung for a tight fit
- Non-coring needle
- With Luer lock connector
- Ergonomic handle
- 100 cm aspiration line




Order No.	Description	Needle Length (cm)	Needle Gauge	Qty. / Box
A-OVM-SL16/35	ACE Single Lumen Ovum Pickup Needle	35	16	10
A-OVM-SL17/35	ACE Single Lumen Ovum Pickup Needle	35	17	10
A-OVM-SL18/35	ACE Single Lumen Ovum Pickup Needle	35	18	10
*All products may not be availa	able in all countries. Please contact our sales department for details			








ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line



The ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line is used for laparoscopic or ultrasoundguided transvaginal retrieval of oocytes from ovarian follicles.

- Highly echogenic tip
- Silicone bung for a tight fit
- Non-coring needle
- With Luer lock connector
- Ergonomic handle
- 100 cm aspiration line
- Alternate aspiration and flushing made possible through a single channel
- Multiple operation functionality by using stopcock valve



Description	Needle Length (cm)	Needle Gauge	Qty. / Box
ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line	35	16	10
ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line	35	17	10
ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line	35	18	10
	ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line	Description Length (cm) ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line 35 ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line 35	DescriptionLength (cm)GaugeACE-FL Single Lumen Ovum Pickup Needle with Flushing Line3516ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line3517







ACE-M Manual Single Lumen Ovum Pickup Needle



The ACE-M Manual Single Lumen Ovum Pickup Needle is used for laparoscopic or ultrasound guided transvaginal retrieval of oocytes from follicles.

- When used with ultrasonic equipment, echo tip enhances visualization
- Handle designed for greater control during ovum retrieval
- Luer lock connector for connecting a syringe
- The amount of suction applied is manually controlled

Tri-facet geometry with highly echogenic tip



Order No.	Description	Needle Length (cm)	Needle Gauge	Qty. / Box
A-OVM-SL16/35-M	ACE-M Manual Single Lumen Ovum Pickup Needle	35	16	10
A-OVM-SL17/35-M	ACE-M Manual Single Lumen Ovum Pickup Needle	35	17	10
A-OVM-SL18/35-M	ACE-M Manual Single Lumen Ovum Pickup Needle	35	18	10
*All products may not be availab	le in all countries. Please contact our sales denartment for details			

*All products may not be available in all countries. Please contact our sales department for details.






WOMEN'S HEALTH




OVUMPICK Double Lumen Ovum Pickup Needle



The OVUMPICK Double Lumen Ovum Pickup Needle is used for laparoscopic or ultrasound-guided transvaginal retrieval of oocytes from ovarian follicles.

- Highly echogenic tip
- Silicone bung for a tight fit
- Non-coring needle
- With Luer lock connector
- The outer lumen allows flushing of follicles




Order No.	Description	Needle Length (cm)	Needle Gauge	Qty. / Box
OVM-DL16/35	OVUMPICK Double Lumen Ovum Pickup Needle	35	16	10
OVM-DL17/35	OVUMPICK Double Lumen Ovum Pickup Needle	35	17	10
*All products may not be availa	able in all countries. Places contact our cales department for details			







EMBRYO TRANS Embryo Transfer Catheter



The EMBRYO TRANS Embryo Transfer Catheter is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- A soft inner catheter with smooth tip
- Marking on inner catheter and outer sheath
- Shapeable outer sheath
- Stylet for outer sheath available separately



Stylet for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	

Order Ne	Description	Catheter
Order No.	Description	Length
EB-ETC23	EMBRYO TRANS Embryo Transfer Catheter	23
EB-ETC18	EMBRYO TRANS Embryo Transfer Catheter	18

STY23	Stylet for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	23	10
STY18	Stylet for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	18	10
STY23-OS	Stylet with outer sheath for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	23	10
STY18-OS	Stylet with outer sheath for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	18	10
*All products may n	ot be available in all countries. Please contact our sales department for details.	/	



WOMEN'S HEALTH

Qty. /


Box

10

10



ECHO TRANS Embryo Transfer Catheter



The ECHO TRANS Embryo Transfer Catheter is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- Echogenic inner catheter facilities enhanced visibility under ultrasound
- A soft inner catheter with a smooth tip
- Marking on inner catheter and outer sheath
- Shapeable outer sheath
- Stylet for outer sheath available separately



Stylet for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter

Order No.	Order No. Description			
ET-ETC23	ECHO TRANS Embryo Transfer Catheter	23	10	
ET-ETC18	ECHO TRANS Embryo Transfer Catheter	18	10	
STY23	Stylet for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	23	10	
STY18	Stylet for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	18	10	
STY23-OS	Stylet with outer sheath for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	23	10	
STY18-OS	Stylet with outer sheath for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	18	10	
*All products may po	the available in all countries. Please contact our sales department for details			

*All products may not be available in all countries. Please contact our sales department for details.






WOMEN'S HEALTH

www.allwinmedical.com




TWINKLE EMBRYO TRANS Embryo Transfer Catheter



The TWINKLE EMBRYO TRANS Embryo Transfer Catheter is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- A soft inner catheter with fused echogenic band
- Marking on inner catheter and outer sheath
- Shapeable outer sheath
- Stylet for outer sheath available separately



Stylet for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter

Order No.	Description	Catheter Length	Qty. / Box
TEB-ETC23	TWINKLE EMBRYO TRANS Embryo Transfer Catheter	23	10
TEB-ETC18	TWINKLE EMBRYO TRANS Embryo Transfer Catheter	18	10
STY23	Stylet for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	23	10
STY18	Stylet for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	18	10
STY23-OS	Stylet with outer sheath for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	23	10
STY18-OS	Stylet with outer sheath for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	18	10
*All products may no	ot be available in all countries. Please contact our sales department for details.		





Trial Transfer Catheter



The Trial Transfer catheter is a closed-ended catheter that is only to be used to determine the passage through the cervix prior to embryo transfer accessible to Allwin's Embryo Transfer Catheter.

- Allows for confident positioning of the outer sheath before introducing the embryo loaded catheter
- Same soft material and smooth tip as catheter to minimize patient trauma
- Centimetre depth marking for confident and precise catheter placement
- Limulus Amebocyte Lysate (LAL) tested to ensure the product is non-pyrogenic




Stylet for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter

Order No.	Description	Catheter Length	Qty. / Box
TT-EB-ETC23	Trial Transfer Catheter for Embryo Trans	23	10
TT-EB-ETC18	Trial Transfer Catheter for Embryo Trans	18	10
TT-ET-ETC23	Trial Transfer Catheter for Echo Trans	23	10
TT-ET-ETC18	Trial Transfer Catheter for Echo Trans	18	10
TT-TEB-ETC23	Trial Transfer Catheter for Twinkle Embryo Trans	23	10
TT-TEB-ETC18	Trial Transfer Catheter for Twinkle Embryo Trans	18	10






BULB TRANS ULTRA Embryo Transfer Catheter - Bulb Tip



The BULB TRANS ULTRA Embryo Transfer Catheter - Bulb Tip is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- Firm outer sheath
- A supported soft inner catheter with a smooth tip
- Marking on inner steel pipe and outer sheath
- Bulb shaped tip of the outer sheath
- Silicone depth positioner for precise placement
- The bulb-shaped angled outer sheath eases passage through the cervix




Stylet for BULB TRANS Embryo Transfer Catheter-Bulb Tip

Description		Qty. / Box
BULB TRANS ULTRA Embryo Transfer Catheter - Bulb Tip		10
Stylet for Bulb Trans Ultra / Bulb Trans Star Embryo Transfer Catheter - Bulb Tip		10
	BULB TRANS ULTRA Embryo Transfer Catheter - Bulb Tip	BULB TRANS ULTRA Embryo Transfer Catheter - Bulb Tip






BULB TRANS STAR Embryo Transfer Catheter - Bulb Tip



The BULB TRANS STAR Embryo Transfer Catheter - Bulb Tip is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- Firm outer sheath
- Soft inner catheter
- Marking on inner steel pipe and outer sheath
- Bulb shaped angled tip of the outer sheath
- Silicone depth positioner for precise placement
- The catheter is curved to facilitate insertion
- Echogenic band on the inner catheter to improve the visibility under ultrasound
- Echogenic tip of the inner catheter to improve the visibility under ultrasound




Stylet for BULB TRANS Embryo Transfer Catheter-Bulb Tip

Description	Qty. / Box
BULB TRANS STAR Embryo Transfer Catheter - Bulb Tip	10
Stylet for Bulb Trans Ultra / Bulb Trans Star Embryo Transfer Catheter - Bulb Tip	10
	BULB TRANS STAR Embryo Transfer Catheter - Bulb Tip





PRO EMBRYO TRANS Embryo Transfer Catheter



Soft Obturator for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter

<u></u>		
Stylet	for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter	
Order No.	Description	Qty. / Bo
PEB-ETC	PRO EMBRYO TRANS Embryo Transfer Catheter	10
PEB-ETC-SB	PRO EMBRYO TRANS Embryo Transfer Catheter with Soft Obturator	10
PEB-ETC-STY	PRO EMBRYO TRANS Embryo Transfer Catheter with Stylet	10
SB-PE	Soft Obturator for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter	10

Stylet for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter

*All products may not be available in all countries. Please contact our sales department for details.




STY-PE

10



PRO ECHO TRANS Embryo Transfer Catheter



Soft Obturator for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter

Stylet for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter


Order No.	Description	Qty. / Box
PET-ETC	PRO ECHO TRANS Embryo Transfer Catheter	10
PET-ETC-SB	PRO ECHO TRANS Embryo Transfer Catheter with Soft Obturator	10
PET-ETC-STY	PRO ECHO TRANS Embryo Transfer Catheter with Stylet	10
SB-PE	Soft Obturator for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter	10
STY-PE	Stylet for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter	10








TWINKLE PRO EMBRYO TRANS Embryo Transfer Catheter





Soft Obturator for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter

	-	-	-	-	di l
_					3
	 -	-	-	-	- N


Stylet for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter

Description	Qty. / Box
TWINKLE PRO EMBRYO TRANS Embryo Transfer Catheter	10
TWINKLE PRO EMBRYO TRANS Embryo Transfer Catheter with Soft Obturator	10
TWINKLE PRO EMBRYO TRANS Embryo Transfer Catheter with Stylet	10
Soft Obturator for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter	10
Stylet for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter	10
	TWINKLE PRO EMBRYO TRANS Embryo Transfer Catheter TWINKLE PRO EMBRYO TRANS Embryo Transfer Catheter with Soft Obturator TWINKLE PRO EMBRYO TRANS Embryo Transfer Catheter with Stylet Soft Obturator for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter






INSEMINA Intra Uterine Insemination Catheter




The INSEMINA Intra Uterine Insemination Catheter (Open Tip & Close Tip) is intended for the introduction of washed spermatozoa into the uterine cavity.

- A soft uterine catheter with firm support
- Integrally molded hub
- Silicone depth positioner for precise placement
- Open tip as well as closed tip (with two lateral holes) catheters available



Order No.	Description	Qty. / Box
I-IUI	INSEMINA Intra Uterine Insemination Catheter	25
I-IUI-O	INSEMINA Intra Uterine Insemination Catheter (Open Tip)	25
*All products may not be	available in all countries. Please contact our sales department for details	











IUI CURVED Intra Uterine Insemination Catheter Curved



The IUI CURVED Intra Uterine Insemination Catheter Curved (Open Tip & Close Tip) is intended for the introduction of washed spermatozoa into the uterine cavity.

- A soft uterine catheter with firm support
- Integrally molded hub
- Silicone depth positioner for precise placement
- Curve shaped catheter
- Open tip as well as closed tip (with two lateral holes) catheters available



Order No.	Description	Qty. / Box
I-IUI-C	CURVED Intra Uterine Insemination Catheter	25
I-IUI-C-O	CURVED Intra Uterine Insemination Catheter (Open Tip)	25
\frown	\frown	\sim
HSSA		


Close Tip



WOMEN'S HEALTH




TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet



The TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet (Open Tip & Close Tip) is intended for the introduction of washed spermatozoa into the uterine cavity.

- A soft uterine catheter with firm support
- Integrally molded hub
- Silicone depth positioner for precise placement
- Malleable stylet to aid catheter insertion
- Open tip as well as closed tip (with two lateral holes) catheters available




Order No.	Description	Qty. / Box
T-IUI	TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet	25
T-IUI-O	TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet (Open Tip)	25
*All products may not be	available in all countries. Please contact our sales department for details.	








FLEXI Intra Uterine Insemination Catheter



The FLEXI Intra Uterine Insemination Catheter (Open Tip & Close Tip) is intended for the introduction of washed spermatozoa into the uterine cavity.

- Shapeable outer sheath for easy insertion
- Centimetre marking on a catheter for precise placement
- Stylet available separately
- Open tip as well as closed tip (with two lateral holes) catheters available



Order No.	Description	Qty. / Box
F-IUI	FLEXI Intra Uterine Insemination Catheter	25
F-IUI-O	FLEXI Intra Uterine Insemination Catheter (Open Tip)	25
STY18	Stylet for FLEXI IUI Intra Uterine Insemination Catheter	25







STUMPY Intra Uterine Insemination Cannula



The STUMPY Intra Uterine Insemination Cannula is intended for the introduction of washed spermatozoa into the uterine cavity.

- Rounded tip
- Luer lock
- Flexible cannula
- Small size catheter




Order No.	Description	Qty. / Box
ST-IUI	STUMPY Intra Uterine Insemination Cannula	25
*All products may not	be available in all countries. Please contact our sales department for details.	






SALPINX CATH Fallopian Tube Recanalization Catheter



The SALPINX CATH Fallopian Tube Recanalization Catheter is used for fallopian tube recanalization under radiological fluoroscopic control.


- Preshaped 5 Fr catheter
- · Centimetre marking on the catheter for precise placement
- Malleable stylet
- Hydrophilic guidewire available separately




Description	Qty. / Box
SALPINX CATH Fallopian Tube Recanalization Catheter with Malleable Stylet	5
SALPINX CATH Fallopian Tube Recanalization Catheter with Malleable Stylet and Hydrophilic Guidewire	1
-	SALPINX CATH Fallopian Tube Recanalization Catheter with Malleable Stylet SALPINX CATH Fallopian Tube Recanalization Catheter




HYSTERO CATH Hysterosalpingography Catheter





The HYSTERO CATH Hysterosalpingography Catheter is used for injection of contrast media into the uterine cavity and fallopian tube(s) to evaluate the patency of the fallopian tubes.

- Strong low profile 7 Fr balloon
- Closed tip with side hole and open tip HSG catheter available
- Stopcock designed for locking the balloon channel and regulating the flow of contrast media.
- Stylet provided for added rigidity during the introduction
- 3 ml syringe included



Description	Qty. / Box
HYSTERO CATH Hysterosalpingography Catheter	5
HYSTERO CATH Hysterosalpingography Catheter Open Taper Tip	5
	HYSTERO CATH Hysterosalpingography Catheter


*All products may not be available in all countries. Please contact our sales department for details.

Caution:

This product contains natural rubber latex which may cause allergic reaction in some individuals.



THE SAMPLER Endometrial Biopsy Catheter



THE SAMPLER Endometrial Biopsy Catheter is used for performing histological tissue biopsy and extracting samples of uterine menstrual content. Used for obtaining endometrial tissue.

- Centimetre marking provided on the catheter for precise placement
- Smooth rounded tip

Smooth rounded tip




Order No.	Description	Qty. / Box
SEBC	THE SAMPLER Endometrial Biopsy Catheter	25
* 4 11	envitable in all according. Discourse and an environment for datable	

*All products may not be available in all countries. Please contact our sales department for details.

PRO-VAC Endometrial Biopsy Sampler The PRO-VAC Endometrial Biopsy Sampler is used for performing histological tissue biopsy and extracting samples of uterine menstrual content. Used for obtaining endometrial tissue. · Centimetre marking provided on the catheter for precise placement • Smooth rounded curved tip Locking syringe Smooth rounded curved tip Centimeter marking provided Locking syringe Order No. Description Qty. / Box **PVES PRO-VAC Endometrial Biopsy Sampler** 25 *All products may not be available in all countries. Please contact our sales department for details.




SURESTRIP Vitrification Straws (Open System Vitrification)



The SURESTRIP Vitrification Straw is a cryopreservation storage device. It is used in the vitrification procedure to contain and maintain oocytes and embryos.

- High survival and development rates
- Secure coupling
- Marking on tip extreme for easy visibility under the microscope
- Adequate flat surface for writing
- Tip width ensuring placement of minimum solution
- 1 cell MEA \ge 80% expand blastocysts within 96h
- Endotoxin LAL tested
- Sterility: SAL 10 $^{\circ}$ Gamma sterilized
- Vitrification of oocytes and embryos in all stages of development
- Available in 5 colours



Colour	White	Red	lvory	Blue	Yellow	Qty. / Box
Order No.	SUS-W	SUS-R	SUS-I	SUS-B	SUS-Y	10





SURELOCK Vitrification Straws



The SURELOCK Vitrification Straw is a cryopreservation storage device. It is used in the vitrification procedure to contain and maintain oocytes and embryos.

- Concave tip for exact and secure placement of specimen under microscope
- Secure air tight coupling
- Semi-transparent body and cap
- The arrow on top surface ensuring orientation during use
- Easy visibility under microscope
- 1 cell MEA \geq 80% expand blastocysts within 96h
- Endotoxin LAL tested
- + Sterility: SAL 10 $^{\circ}$ Gamma sterilized
- Vitrification of 1-cell stage of embryos
- Available in 5 colours



specimen under microscope

Colour	Pink	Yellow	Green	Blue	Violet	Qty. / Box
Order No.	SUL-P	SUL-Y	SUL-G	SUL-B	SUL-V	10


*All products may not be available in all countries. Please contact our sales department for details.



WOMEN'S HEALTH




SUREHOOK Vitrification Straws



The SUREHOOK Vitrification Straw is a cryopreservation storage device. It is used in the vitrification procedure to contain and maintain oocytes and embryos.

- Hook shaped tip design
- Secure coupling
- Adequate flat surface for writing
- Hook width ensuring placement of minimum solution
- Flexible tip
- 1 cell $\dot{MEA} \ge 80\%$ expand blastocysts within 96h
- Endotoxin LAL tested
- Sterility: SAL 10⁻⁶ Gamma sterilized
- Vitrification of oocytes and embryos
- Available in 5 colours



Colour	White	Red	Green	Blue	Yellow	Qty. / Box
Order No.	SUH-W	SUH-R	SUH-G	SUH-B	SUH-Y	10




allwin[®] Medical Devices

ASPI Aspiration Needle




The ASPI Aspiration Needle is used for laparoscopic or ultrasound-guided transvaginal aspiration.

- Bevel tip
- Centimetre marking provided for precise placement
- Transparent luer lock connector



Order No.	Description	Size cm	Gauge	Qty. / Box
ASPI 18/30	ASPI Aspiration Needle	30	18	10
ASPI 22/30	ASPI Aspiration Needle	30	22	10





Female Reproductive System

Balance Medical Devices

US Sales

sales.usa@allwinmedical.com
► +1 714-572-1709

India Sales

sales.india@allwinmedical.com
► +91 99870-84563

Worldwide Sales (Excluding India/US)

≤ sales@allwinmedical.com

\$ +91 95949-81581