

SURELOCK

Vitrification Straws



SURELOCK

Vitrification Straws

The SURELOCK Vitrification Straw is a cryopreservation storage device. It is used in the vitrification procedure to contain and maintain oocytes and embryos.



Marking for easy visualization during insertion



Concave tip for exact and secure placement of specimen under microscope



Flat surface wide body for writing

Colour	Pink	Yellow	Green	Blue	Violet	Qty. / Box
Order No.	SUL-P	SUL-Y	SUL-G	SUL-B	SUL-V	10

*All products may not be available in all countries. Please contact our sales department for details.

US Sales

✉ sales.usa@allwinmedical.com
☎ +1 714-572-1709

India Sales

✉ sales.india@allwinmedical.com
☎ +91 99870-84563

Worldwide Sales (Excluding India/US)

✉ sales@allwinmedical.com
☎ +91 95949-81581