

TRACKABLE IUI



Intra Uterine Insemination Catheter

allwin[®]
Medical Devices
www.allwinmedical.com

TRACKABLE IUI

Intra Uterine Insemination Catheter with Malleable Stylet


The TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet (Open Tip & Close Tip) is intended for the introduction of washed spermatozoa into the uterine cavity.



Silicone depth positioner



Integrally molded hub



Firm Support



With malleable stylet

Order No.	Description	Qty. / Box
T-IUI	TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet	25
T-IUI-O	TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet (Open Tip)	25

*All products may not be available in all countries. Please contact our sales department for details.



Close Tip



Open Tip

US Sales

✉ sales.usa@allwinmedical.com
☎ +1 714-572-1709

India Sales

✉ sales.india@allwinmedical.com
☎ +91 99870-84563

Worldwide Sales (Excluding India/US)

✉ sales@allwinmedical.com
☎ +91 95949-81581