

# TWINKLE EMBRYO TRANS



Embryo Transfer Catheter

# TWINKLE EMBRYO TRANS

## Embryo Transfer Catheter


The TWINKLE EMBRYO TRANS Embryo Transfer Catheter is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.




Stylet for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter



Fused  
echogenic Band



Marking on  
inner catheter



Orientation  
arrow on hub

Order No.	Description	Catheter Length	Qty. / Box
TEB-ETC23	TWINKLE EMBRYO TRANS Embryo Transfer Catheter	23	10
TEB-ETC18	TWINKLE EMBRYO TRANS Embryo Transfer Catheter	18	10
STY23	Stylet for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	23	10
STY18	Stylet for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	18	10
STY23-OS	Stylet with outer sheath for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	23	10
STY18-OS	Stylet with outer sheath for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	18	10

\*All products may not be available in all countries. Please contact our sales department for details.



### US Sales

✉ sales.usa@allwinmedical.com  
☎ +1 714-572-1709

### India Sales

✉ sales.india@allwinmedical.com  
☎ +91 99870-84563

### Worldwide Sales (Excluding India/US)

✉ sales@allwinmedical.com  
☎ +91 95949-81581