

TWINKLE PRO EMBRYO TRANS



Embryo Transfer Catheter

TWINKLE PRO EMBRYO TRANS

Embryo Transfer Catheter

The TWINKLE PRO EMBRYO TRANS Embryo Transfer Catheter is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

TWINKLE PRO EMBRYO TRANS Embryo Transfer Catheter is also available with compatible stylet or soft obturator which assist the guide catheter in gaining uterine access in the embryo transfer procedure.




Soft Obturator for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter



Stylet for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter



Soft inner catheter with fused echogenic band



Silicone depth positioner



Centimetre marking



Orientation arrow on hub

| Order No. | Description | Qty. / Box |
|-------------|--|------------|
| TPEB-ETC | PRO EMBRYO TRANS Embryo Transfer Catheter | 10 |
| TPEB-ETC-SB | PRO EMBRYO TRANS Embryo Transfer Catheter with Soft Obturator | 10 |
| TEB-ETC-STY | PRO EMBRYO TRANS Embryo Transfer Catheter with Stylet | 10 |
| SB-PE | Soft Obturator for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter | 10 |
| STY-PE | Stylet for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter | 10 |

*All products may not be available in all countries. Please contact our sales department for details.



US Sales

✉ sales.usa@allwinmedical.com
☎ +1 714-572-1709

India Sales

✉ sales.india@allwinmedical.com
☎ +91 99870-84563

Worldwide Sales (Excluding India/US)

✉ sales@allwinmedical.com
☎ +91 95949-81581